



Compliance. Safety. Function.

FloorBed¹
by Accora

FloorBed[®] 1

by Accora



Full height



Comfort position



Care position

FLOORBED[®] 1 TECHNICAL SPECIFICATION

Dimensions:	36.4" x 90.4"	Lifting range:	2.8" – 25.6" (excl. mattress)
Mattress size:	Standard configuration: 36.0" x 80.0" Extended configuration*: 36.0" x 84.0"	Profiling functions:	Backrest with Auto-regression, Kneebreak, Anti-Trendelenburg, Trendelenburg
Safe working load:	408 lbs	Part number:	NSB-0-FL1-200US
Max patient weight:	330 lbs		

*Mattress platform extension infill piece MUST be used if the mattress platform is extended.

COLOR OPTIONS



CHERRY



OAK

FloorBed[®] 1Plus

by Accora

Bariatric features

- 42" Mattress Platform
- 573 lbs Safe Working Load



FLOORBED[®] 1 PLUS TECHNICAL SPECIFICATION

Dimensions:	42.3" x 90.4"	Lifting range:	3.1" – 25.6" (excl. mattress)
Mattress size:	Standard configuration: 42.0" x 80.0" Extended configuration*: 42.0" x 84.0"	Profiling functions:	Backrest, Kneebreak, Anti-Trendelenburg, Trendelenburg
Safe working load:	573 lbs	Part number:	NSB-0-FL5-100US
Max patient weight:	496 lbs		

*Mattress platform extension infill piece MUST be used if the mattress platform is extended.

Compliance.

Certified compliance with the latest safety standards, including
 ANSI/AAMI ES60601-1
 IEC 60601-1-2
 IEC 60601-2-52



Safety.



Function.



FloorBed® 1 Features



Width adjustable mattress guides



No trailing wires




Full nursing height



Concealed wheels



Backrest and kneebreak functions



Mattress platform height of just 2.8 inches

FloorBed® 1 Accessories



Head and foot board bumpers



Mattress platform extension – increasing length to 84 inches



Assist bar to assist with bed mobility



Safety mat provides same height as a FloorBed® with 6 inch foam mattress



Mattress pump bracket

Accora

9210 Corporate Blvd. Suite 120, Rockville, MD 20850

T: 301 560 2400

E: information@accora.care

W: accora.care

* and ™ are trademarks of Accora Limited.
As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

